

vacuflex[®]2000

The lymph is a fluid located in the intercellular space, which consists of the fluid being secreted in the arterious side of capillaries and the waste material created during the metabolism of cells. Thus the lymph is in direct contact with the cells and the blood vessels. Any chemical or biological change of the organism will affect lymphatic production and the composition of the intercellular fluid. The lymph is collected in the lymphatic vessels, and filtered by the lymphatic glands located along the lymph routes as a filter system, thus cleaning it of harmful substances.

Improper nutrition, excess uptake of fat and carbohydrates, smoking, or drinking coffee, and use of oral contraceptives influence the work of the lymphatic system - due to the growth of fatty cells, damage to lymphatic vessels, and depositing of waste materials. As a result, many ladies may expect the appearance of the cellulite problem and striae, which represent marked aesthetic imperfection.

Vacuum massage

After analyzing these pathological processes of the cellulite problem, the doctors, natural therapists and engineers involved in developing our treatment devices, developed a treatment method that helps eliminate not only mild, but even more severe forms of the cellulite problem. The beneficial full-bloodedness caused by vacuum enhances metabolism and improves the removal of waste materials.

Owing to the tensile effect of vacuum, the membrane of many fatty cells becomes permeable, the fat flows out through it, and gets into the lymph flow. The vacuum activates the cells of the connective tissue, which produce flexible fibres, thus the skin also becomes more stretched. As a result of the special massaging movements carried out with the help of the vacuum, the speed of lymph flow may be multiplied by 50 even. Thus the arising waste materials (metabolites / catabolites) may be removed through the most natural route possible.

With the help of vacuum you may achieve a complete refreshment of the body, which goes beyond the result of a simple massage treatment. Applying carefully selected active ingredients, you can further increase this effect, and using it in combination with other methods (e.g. ultrasound, interference current or iontophoresis), you can determine a complete body care program.



The reflex therapy and stimulation with aromatic oils

The beneficial effects of vacuum are well known for thousands years, according to the finds from Mesopotamia, the people of that time already used it in 3000 B.C.

The glass accessory placed on the skin - if the air is sucked out from them - can influence the functions of the internal organs through the remote effect of the meridian and reflex points by the vacuum effect. The stimuli may alter the organism through the vegetative regulation, and may increase its adaptation capability.

With the help of the automatic massage programs of the vacuflex® device - using 2, 4 or 6 glass accessories - we can stimulate simultaneously the reflex points.

Beyond these, we can also carry out a massage on the palm or the sole with the handpieces. By placing the glass tubes on the proper area we can achieve a pleasant muscle relaxing in order to relax the tensed, inflexible muscles.

The treatments could be combined with the stimulating and soothing aromatic oils, with which the different therapies could be made even more efficient. The cupping accessories can be placed on the back after a massage with essential oils.

We can apply prepared, made depending on the individual characteristics or freshly made - with strengthening, soothing or drainage effect - combinations of oils.

The effects of vacuum:

- local blood circulation improves
- livelier metabolism
- relief of muscle spasms
- oxygen intake of cells increase
- stimulation of the immune system
- mobilization of fatty deposits
- analgesic (pain relieving) effect
- improved lymph flow



Easy contact

The comfort offered by vacuflex® 2000 becomes obvious already when you connect the handpiece to it. The treating tube is fitted with the new easy contact socket, and there is a hermetic connection immediately. So you can be sure that the vacuum produced by the device will be utilized completely, and the treating tubes wouldn't be damaged even by constant use.

vacucontrol®

The high performance compressor of the vacuflex® 2000 device ensures the quick attainment of the preset sucking power, thus we can carry out the treatment with precise movements, without any unnecessary waiting. The vacuflex® 2000 device continuously monitors the sucking force during the treatment, and adjusts it according to the output selected on the front panel. The continuously controlled output of the vacuum ensures that you will always work with a vacuum suitable for the given part of the body. The vacucontrol® guarantees a successful treatment free of inconveniences, so that your client will be really content.



Without any trouble

If you start working with the vacuflex® 2000 device, you will feel immediately: this instrument has been designed for you. You don't need to trouble with inflexible tubes and unwieldy hadpieces: the light and easy to handle handpieces, the transparent accessories made from tempered glass, the graceful connecting tubes radiate comfort and elegance in all treating position.

Our designer has seen the vacuflex® 2000 in dreams for you, now the success depends only on you. Without any trouble.

Technical details

Output parameter of vacuum:

20-400 mBar in 20 mBar steps

Dimensions: 255 x 245 x 140 mm

Weight: 3 kg



Basic equipment beyond your wildest expectations

The vacuflex® 2000 deserves a distinct attention not only for its technical layout, but also for the services provided by it.

The vacuum massage and cellulite therapy set, which are part of the basic equipment, include a handpiece with 3 glass accessories of different size, and a fast connection tube for the easy contact socket. Furthermore, the set contains 2 handpieces and a distributor for face treatment carried out with two hands simultaneously.

In addition to vacuum massage treatment of the face and body, there are also various reflex therapy programs, which serve a full-scope treatment concept with the help of the optionally available reflex therapy set.

