

Ultrapeel® Pepita System

Easy to use

Equipped with No clogging system

Sterile Disposables accessories

Patented 45° V shaped light plastic tips

Patented Disposable sealed waste bottle

Patented Disposable sealed and sterile clean crystals bottle

Power adjustment features

External compressor connection

Suitable for:

-Facial and body peeling

-Aging skins

-Dry, oily, dull, congested and blemished skins

-Open pores

-Wrinkles

-Acne

-Scars

-Hyper pigmentation and Age spots

-Stretch Marks reduction

Pre-treatment for further applications (Laser and Chemical Peel)



The Ultrapeel Pepita System is an innovative approach for skin resurfacing.

The practitioner can carry out a wide range of aesthetic and para medical procedures, performing non invasive and progressive treatments. Most patients can continue with a normal routine while undergoing a course of treatments.

The equipment uses a controlled flow of inert sterile micro crystals, which are vacuumed across the skin removing dead and damaged skin cells. The exfoliation of the skin comes through a jet of sterile inert corundum crystals, which interact with the tissues through a patented 45 deg. V-shaped sterile disposable tip, in a homogeneous and efficient way.

The treatment is carried out in a safe and efficient way, satisfying any sort of customer's expectations.

The Ultrapeel® Pepita uses a complete sterilized and disposable delivery kit in order to avoid any kind of contamination.

It means that it is the only Microdermabrasion System, together with the Aesthi Peel® and the other Ultrapeel® Systems, granting sterile operations.

The high reliability of this equipment is based on its exclusive characteristics.

It is in fact equipped with an anti clogging system, and the high durability is also guaranteed by top quality components and its exclusive filtering capability. The light easy to use accessories and the power regulation make the Ultrapeel Pepita to be the standard in the microdermabrasion field. It was also cleared by the FDA on December 1996, together with the Ultrapeel® Crystal System.

ACCESSORIES AND DISPOSABLES



*Pre-Packed, Pre-Sterilized,
disposable delivery Kit,
accessories and corundum
#30000285 Corundum
bottle set
#30010015 Disposable
Tubing kit*



*45 Deg. V-Shaped
Patented disposable
abrading tip
#10010142*

TECHNICAL SPECIFICATIONS

POWER SUPPLY:

Supply Voltage: 110 and 220 V,
single phase

Frequency : 50 / 60 Hz

Power absorption : 300 VA

Fuses: 2x3,15 AT at 110V, 60Hz
2x1,6 AT at 220V, 50Hz

PERFORMANCES:

Suction power: from 0.2 to 0.7
bar (at sea level)

Abrading media:

High purity medical grade
corundum powder

PHYSICAL:

Dimensions: 30x30x45cm.
12x12x18 in

Weight: 18 kg / 39.lb

STERILIZATION

Sterilization is of paramount importance when the operator is dealing with different patients and skin conditions.

Cross Contaminations and debris residua may affect the patient if the Microdermabrasion system is not perfectly clean and sterile.

Sterilized disposables make sure that nothing happens to the patient and operator.

Only Mattioli Engineering guarantees safe and sterile operations through its disposable and consumables accessories

Mattioli Engineering - A new science in Skin Care®

Mattioli Engineering Italia SPA

Villa dell'Ombrellino
Piazza di Bellosguardo 11
50124 FIRENZE—ITALY
Phone.: +39 055 221603
Fax: +39 055 221735
www.matteng.com
info@mattioliengineering.it

**Mattioli Engineering Corporation**

8300 Greensboro Drive
Suite 800—McLean VA 22102
Toll Free: 877 628 8364
Phone (703) 312 6000
Fax: (703) 243 9139
www.matteng.com
info@matteng.com