

Ultrapeel® II plus — Microdermabrasion Unit

Easy to use and no clogging

New Features:

Digital Crystal flow adjustment

Digital Pressure adjustment

Digital Vacuum adjustment

Inbuilt Compressor

Sterile Disposables accessories

Patented 45° V shaped light handpiece

Patented Disposable sealed waste bottle

Patented Disposable sealed and sterile clean crystals bottle

Transderm® Compatible
Ultraprobe™ Compatible

Suitable for:

- Facial and body peeling*
- Ageing skins*
- Dry, oily, dull, congested and blemished skins*
- Open pores & Wrinkles*
- Acne*
- Scars*
- Hyper pigmentation and Actinic chertosis*
- Stretch Marks reduction*
- Tattoo removal*
- General dermabrasion.*

Pre-treatment for further applications (Laser and Chemical Peel)



The latest born into the Ultrapeel® Family of products.

Ultrapeel® II **plus** microdermabrasion system represents the most innovative and technologically advanced medical Micordermabrasion system suitable for dermatologists, Plastic and Cosmetic Surgeons, Medical doctors. Ultrapeel® II **plus** allows users to perform a huge range of treatments in a complete safe , Hygienic ,sterile environment with corundum crystals saving, avoiding any cross contamination risks completely and to satisfy customer's expectations of any sort.

Save crystal for Real....

Ultrapeel® II **plus** disposable delivery kit exfoliates the skin through a jet of inert corundum crystals through the patented 45 deg. V-shaped handpiece in a uniform and effective way , abrading the skin in a gentle but efficient way.

No Clogging and durable, Ultrapeel® II **plus** gives users the opportunity to **save crystals** by the latest **Crystal flow regulation** feature on board. In addition to this the **maximum vacuum power setting** and adjustment feature makes the machine able to be configured for any type of use in complete safety for either patients and users taking care of the user and the crystal consumption .

In built compressor for critical applications....

The foot-pedal that comes up with Ultrapeel® II **plus** , enables air compression for deeper treatments bombarding the patient's skin with sterile corundum and then vacuum it out leaving the skin free from any residua.

The most reliable and versatile machine on the market

The light easy to use accessories, the sterilization features and the power regulation make the Ultrapeel® II **plus** to set **standards of care** in the Medical field for any patient and any skin condition. Ultrapeel® have been the very first microdermabrasion systems to be cleared, by the FDA on December 1996.

Ultrapeel® II **plus** is compatible for the **Transderm® Meso** Transdermal delivery system the very first FDA cleared alternative to injection. (see the *Transderm lonto data sheet for details*).

Ultrapeel® II **plus** is already compatible for the new ULTRAPROBE™ Universal pre Filled , pre sterilized, disposable handpiece (see the *Ultraprobe™ data sheet for details*).

ACCESSORIES AND DISPOSABLES



*Pre-Packed, Pre-Sterilized,
disposable delivery Kit,
accessories and corundum
#30000385 Corundum
bottle set
#30010010 Disposable
Tubing kit*



*45 Deg. V-Shaped
Patented disposable
abrading tip
#10010142*

**Mattioli Engineering Italia SPA**

P.zza di Bellosguardo 11
50124 FIRENZE—ITALY
Phone.: +39 055 221603
Fax: +39 055 221735
info@mattioli-engineering.it

Mattioli Engineering Corporation

830 Greensboro Drive
Suite 800—McLean VA 22102
Toll Free: 877 628 8364
Phone (703) 312 6000
Fax: (703) 243 9139
www.matteng.com
info@matteng.com

TECHNICAL SPECIFICATIONS

POWER SUPPLY:

Supply Voltage: 220 ,single phase
Frequency: 50 Hz
Power absorption 440VA
Fuses: 2 x 2 AT at 220V, 50Hz
Class (according to IEC 601-1)1 type
BF for Dermoelectroporation
Class 1 type B for Dermabrasion
Extn Connector: 9Vdc 1A max
Compression: 2,5 bar max

PERFORMANCES:

Suction power: from 0.2 to 0.7 bar
(below 10000 mt altitude)
Abrading media:
High purity medical grade corundum
powder

PHYSICAL:

Dimension : 32x35x55 cm
Weight 20 kg

STERILIZATION AND SAFETY

Sterilization is of paramount importance when the operator is dealing with different Patients and skin conditions.

Cross Contaminations and debris residua may affect the patient if the microdermabrader is not perfectly cleaned and Sterile.

Sterilized disposables make sure that nothing will happen to the patient and to the operator at all.

Only Mattioli Engineering guarantees safe and sterile operations thru its disposable accessories and consumables.

Mattioli Engineering - A new science in Skin Caretm